

Feed Water Requirements

Model	AquaRO ²	EDI ²
type of feed water	potable water / tap water	potable water / tap water
feed water conductivity	< 1400 μ S/cm	< 1400 μ S/cm
inlet pressure	1.5 to 6 bar	< 5 bar
free chlorine	< 0.1 mg/L	< 0.02
Silt Density Index (SDI)	< 3	< 3
pH	3 to 9	3 to 9
temperature	5 to 25 °C	5 to 25 °C
iron manganese	n. a. n. a.	< 0.01 mg/L < 0.01 mg/L
sulfate	n. a.	< 2 mg/mL