

Technical Specifications



WATER
PURIFIER

Model	AquaRO ²	AquaRO ² DI	AquaRO ² 8 P10 P10 DI	EDI ²
pure water quality	< 10 µS/cm	< 1 µS/cm; Type II	< 10 µS/cm < 1 µS/cm	< 0.2 µS/cm; Type II
resistivity	> 0.1 MΩ·cm	> 1 MΩ·cm	> 0.1 MΩ·cm > 1 MΩ·cm	> 5 MΩ·cm
total organic carbon (TOC)	< 50 ppb	< 40 ppb	< 40 ppb	< 20 ppb
productivity rate	8 L/h to 40 L/h	8 L/h to 20 L/h	8 L/h	10 l/h to 40 L/h
bacteria	< 50 cfu/mL	< 50 cfu/mL	< 50 cfu/mL	< 10 cfu/mL

dimensions, weight, power (except AquaRO² 40 and EDI² 40):

AquaRO² : 504 x 340 x 535 mm, 16 to 20 kg, 110 V - 230 V

EDI² : 504 x 340 x 535 mm, > 40 kg, 110 V - 230 V

Peybord